



Effervescent Cleaning and
Disinfecting Tablets for
Hospital and Institutional Use

SPORICIDAL

ACTIVE INGREDIENT*:

Hypochlorous Acid.....0.074%
Other Ingredients.....99.926%
Total.....100.000%

4306 ppm Available Chlorine

**Percentages may vary slightly based on the active ppm solution.
Refer to the dilution chart for Available Chlorine concentrations.*



**1 Minute
Norovirus**

KEEP OUT OF REACH OF CHILDREN

CAUTION

Follow the Directions for Use on the pesticide label when applying this product. The product in this container is diluted as directed on the pesticide label. This is the ready-to-use dilution of Trust Medical's concentrate Defender.

FIRST AID

IF IN EYES: Hold eye open and rinse slowly and gently with water for 15-20 minutes. If present, remove contact lenses after the first 5 minutes, then continue rinsing eye.

Call a poison control center or doctor for further treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment.

PRECAUTIONARY STATEMENTS TO HUMANS AND DOMESTIC ANIMALS

CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet.

Wear chemical-resistant gloves and safety glasses or face shield when making up solution.

Refill only with this product. Take care to update the expiration date on refillable bottle if concentrate kept past one day and replace bottle if it becomes worn or damaged.

Refer to pesticide label and Safety Data Sheet for additional safety information.

DIRECTIONS FOR USE

General Disinfection Without Preclean:

1. Wear appropriate barrier protection.
2. Prepare a 4306ppm solution (refer to Dilution Chart).
3. Place tablet(s) in cool tap water and let effervesce until dissolved (approx. 3-6 mins).
4. Apply to surface with mop, cloth, wipe etc.
5. Allow surface to remain wet for 4 minutes.
6. Recommend that user prepare a new solution daily. However, solution may be kept for 3 days when in a closed container and dated below. Remaining solution can be poured down drain pipe for drain pipe disinfection.
7. All treated equipment and surfaces that will contact food, feed, or drinking water must be rinsed with potable water before use.

Disinfection For Surfaces Contaminated With C. difficile:

1. Follow steps 1 – 3 above.
2. Apply to pre-cleaned surface with mop, cloth, wipe etc.
3. Allow surface to remain wet for 4 minutes. Allow to air dry
4. Follow steps 6 – 7 above.

REFER TO PESTICIDE LABEL FOR FULL DIRECTIONS FOR USE

| Solution PPM (mg/L) Available Chlorine | Water Volume | 13.1g Tablets | 3.34g Tablets |
|---|-----------------|------------------|------------------|
| 4306 | 32 Ounces | 1 | 4 |
| 4306 | 1 Gallon | 4 | 16 |

**IF SOLUTION IS KEPT PAST ONE DAY
DISCARD & REFILL ON (DAY)**



DISCARD CONTENTS AFTER 3 DAYS

4306 PPM Available Chlorine (C. difficile spores in 4 minutes)



Trust Medical, LLC
1-888-200-0917
www.TrustMedicalCompany.com

EPA Reg. No. 71847-7-92281
EPA Est. No. 71847-IRL-001

Customer Affixed Label